

08/30/2017

Page 1 of 2

Certificate of Analysis

Customer: Gesz Gaal es Sziklas Kft.**Sample Identification:**

Batch #: B-15504

BL ID #: 17-0438

Description: FLAVON PEAK, paste, 00217

Date Received: 08/17/2017

Results:

Analysis	Result	Units
ORAC against peroxy radicals	194	µmole TE/gram
ORAC against hydroxyl radicals	690	µmole TE/gram
ORAC against peroxy nitrite	56	µmole TE/gram
ORAC against super oxide anion	1982	µmole TE/gram
ORAC against singlet oxygen	15,536	µmole TE/gram
ORAC against hypochlorite	368	µmole TE/gram

There are six predominant reactive species found in the body: peroxy radicals, hydroxyl radicals, peroxy nitrite, super oxide anion, singlet oxygen and hypochlorite. ORAC 6.0 provides a measure of the total antioxidant power of a food/nutrition product against the six predominant reactive species.

The ORAC result is expressed as micromole Trolox equivalency (µmole TE) per gram.

References:

- [1] Ou, Boxin, Maureen Hampsch-Woodill, and Ronald L. Prior. "Development and validation of an improved oxygen radical absorbance capacity assay using fluorescein as the fluorescent probe." *Journal of agricultural and food chemistry* 49.10 (2001): 4619-4626.
- [2] Huang, Dejian, et al. "Development and validation of oxygen radical absorbance capacity assay for lipophilic antioxidants using randomly methylated β -cyclodextrin as the solubility enhancer." *Journal of Agricultural and Food Chemistry* 50.7 (2002): 1815-1821.
- [3] Ou, Boxin, et al. "Novel fluorometric assay for hydroxyl radical prevention capacity using fluorescein as the probe." *Journal of Agricultural and Food Chemistry* 50.10 (2002): 2772-2777.
- [4] Dubost, N. Joy, Boxin Ou, and Robert B. Beelman. "Quantification of polyphenols and ergothioneine in cultivated mushrooms and correlation to total antioxidant capacity." *Food Chemistry* 105.2 (2007): 727-735.
- [5] Zhang, Liliang, et al. "Novel high-throughput assay for antioxidant capacity against superoxide anion." *Journal of agricultural and food chemistry* 57.7 (2009): 2661-2667.
- [6] Ou, Boxin, Dejian Huang, and Maureen H. Woodill. "Method for assaying the antioxidant capacity of a sample." U.S. Patent No. 7,132,296. 7 Nov. 2006.

The results shown in this Certificate of Analysis refer only to the sample(s) tested, unless otherwise stated. Attention is drawn to the limitation of liability, indemnification, and jurisdictional issues. This Certificate of Analysis cannot be reproduced, except in full, without prior written permission of Brunswick Laboratories, Inc.. Any misrepresentation, unauthorized alteration, or falsification of the content or appearance of this Certificate of Analysis is unlawful. The Customer agrees and indemnifies Brunswick Laboratories, Inc., its officers, directors, employees, agents, (collectively referred to as Brunswick Labs) and agrees to hold them harmless from any claims, judgments, actions, or expenses of any kind, including attorneys' fees, in the event that Brunswick Labs takes action to correct any such misrepresentation, alteration, forgery, or falsification, understanding that such conduct by Customer or its employees or agents damages the reputation and therefore the business of Brunswick Labs. Such actions may include, without limitation, litigation or announcements to any component of or the entire relevant industry of such misbehavior, as described above.

08/30/2017

Page 2 of 2

Certificate of Analysis

Customer: Gesz Gaal es Sziklas Kft.**Sample Identification:**

Batch #: B-15504

BL ID #: 17-0438

Description: FLAVON PEAK, paste, 00217

Date Received: 08/17/2017

ORAC 6.0 Method Description:

Oxygen Radical Absorbance Capacity (ORAC) tests are among the most acknowledged methods that measure antioxidant scavenging activity against oxygen radicals that are known to be involved in the pathogenesis of aging and many common diseases. ORAC 6.0 consists of six types of ORAC assays that evaluate the antioxidant capacity of a material against six primary reactive oxygen species (ROSs, commonly called "oxygen radicals") found in humans: peroxy radical, hydroxyl radical, superoxide anion, singlet oxygen, peroxynitrite, and hypochlorite. This is a comprehensive panel that evaluates the antioxidant capacity of a material against oxygen radicals.

The ORAC 6.0 tests are based on evaluating the capacity of an interested material to protect a probe (a fluorescent probe or chromagen) from its damage by ROSs. In all ORAC assays, an ROS inducer is introduced to the assay system. The ROS inducer triggers the release of a specific ROS, which would degrade the probe and cause its emission wavelength or intensity change. When an antioxidant material presents in the environment, the antioxidant absorbs the ROS and preserves the probe from degradation. The degree of probe preservation indicates the antioxidant capacity of the material. Trolox is used as the reference standard, and the results are expressed as μ mole Trolox equivalency per gram (or milliliter) of a tested material.

The ORAC 6.0 values obtained from various natural and synthetic materials range widely due to the significantly different natures of these materials including fresh produce, liquid, extract, pure compound, cream, or oil. [Select ORAC 6.0 values of benchmark materials](#) are listed in our website to provide limited scale. Comparisons between materials from similar categories are more informative and valuable.

The results shown in this Certificate of Analysis refer only to the sample(s) tested, unless otherwise stated. Attention is drawn to the limitation of liability, indemnification, and jurisdictional issues. This Certificate of Analysis cannot be reproduced, except in full, without prior written permission of Brunswick Laboratories, Inc.. Any misrepresentation, unauthorized alteration, or falsification of the content or appearance of this Certificate of Analysis is unlawful. The Customer agrees and indemnifies Brunswick Laboratories, Inc., its officers, directors, employees, agents, (collectively referred to as Brunswick Labs) and agrees to hold them harmless from any claims, judgments, actions, or expenses of any kind, including attorneys' fees, in the event that Brunswick Labs takes action to correct any such misrepresentation, alteration, forgery, or falsification, understanding that such conduct by Customer or its employees or agents damages the reputation and therefore the business of Brunswick Labs. Such actions may include, without limitation, litigation or announcements to any component of or the entire relevant industry of such misbehavior, as described above.